

产品名称: KLC1 Rabbit Monoclonal Antibody
产品货号: AMRe84608



产品概述 (Summary)

产品名称 (Production Name)	KLC1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KLC1
别名 (Alternative Names)	Kinesin2; KLC; Klc1; KNS2; KNS2A;;Kinesin light chain 1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q07866.A synthesized peptide derived from human Kinesin light chain 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	65 kDa

研究背景 (Background)


产品名称: KLC1 Rabbit Monoclonal Antibody
产品货号: AMRe84608




Kinesin is a microtubule-associated force-producing protein that may play a role in organelle transport. The light chain may function in coupling of cargo to the heavy chain or in the modulation of its ATPase activity.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of KLC1 expression in HeLa cell lysate.




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: KLC1 Rabbit Monoclonal Antibody
产品货号: AMRe84608



Immunoprecipitate (IP) analysis using the Antibody at 1:50 dilution. (wb at 1:3K dilution)

注意事项 (Note)

For research use only .