

产品名称: APLP2 Rabbit Monoclonal Antibody
产品货号: AMRe84611



产品概述 (Summary)

产品名称 (Production Name)	APLP2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	APLP2
别名 (Alternative Names)	Aplp2; APPH; APPL2; CDEBP;;APLP2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q06481.A synthesized peptide derived from human APLP2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 87 kDa ; Observed MW: 87-140 kDa

研究背景 (Background)

May play a role in the regulation of hemostasis. The soluble form may have inhibitory properties towards coagulation


产品名称: APLP2 Rabbit Monoclonal Antibody
产品货号: AMRe84611




factors. May interact with cellular G-protein signaling pathways. May bind to the DNA 5'-GTCACATG-3'(CDEI box). Inhibits trypsin, chymotrypsin, plasmin, factor XIA and plasma and glandular kallikrein. Modulates the Cu/Zn nitric oxide-catalyzed autodegradation of GPC1 heparan sulfate side chains in fibroblasts.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of APLP2 expression in HeLa cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .