

产品名称: Kallikrein 5 Rabbit Monoclonal Antibody
产品货号: AMRe84621



产品概述 (Summary)

产品名称 (Production Name)	Kallikrein 5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Kallikrein 5
别名 (Alternative Names)	HK5; Kallikrein 5; Kallikrein L2; KLK 5; KLK L2; SCTE;; Kallikrein 5
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9Y337.A synthesized peptide derived from human Kallikrein 5

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 32 kDa ; Observed MW: 37 kDa

研究背景 (Background)

Kallikreins are a subgroup of serine proteases having diverse physiological functions. Growing evidence suggests that many


产品名称: Kallikrein 5 Rabbit Monoclonal Antibody
产品货号: AMRe84621



kallikreins are implicated in carcinogenesis and some have potential as novel cancer and other disease biomarkers. This gene is one of the fifteen kallikrein subfamily members located in a cluster on chromosome 19. Its expression is up-regulated by estrogens and progestins.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Kallikrein 5 expression in MCF-7 cell lysate.

注意事项 (Note)

For research use only .