

产品名称: XPD Rabbit Monoclonal Antibody
产品货号: AMRe84623



产品概述 (Summary)

产品名称 (Production Name)	XPD Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	XPD
别名 (Alternative Names)	BTF2 p80; COFS2; CXPB; EM9; ERCC2; MAG; TFIIH p80; TTD; XPD; XPDC;;XPD
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P18074.A synthesized peptide derived from human XPD

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 87 kDa ; Observed MW: 75 kDa

研究背景 (Background)

ATP-dependent 5'-3' DNA helicase, component of the general transcription and DNA repair factor IIH (TFIIH) core complex,


产品名称: XPD Rabbit Monoclonal Antibody
产品货号: AMRe84623



which is involved in general and transcription-coupled nucleotide excision repair (NER) of damaged DNA and, when complexed to CAK, in RNA transcription by RNA polymerase II.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of XPD expression in A431 cell lysate.

注意事项 (Note)

For research use only .