

产品名称: DMAP1 Rabbit Monoclonal Antibody
产品货号: AMRe84627



产品概述 (Summary)

产品名称 (Production Name)	DMAP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DMAP1
别名 (Alternative Names)	DNA methyltransferase 1 associated protein 1; DNMAP1; DNMTAP1; EAF2; SWC4;; DMAP1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9NPF5.A synthesized peptide derived from human DMAP1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	53 kDa

研究背景 (Background)


产品名称: DMAP1 Rabbit Monoclonal Antibody
产品货号: AMRe84627




Involved in transcription repression and activation. Its interaction with HDAC2 may provide a mechanism for histone deacetylation in heterochromatin following replication of DNA at late firing origins. Can also repress transcription independently of histone deacetylase activity. May specifically potentiate DAXX-mediated repression of glucocorticoid receptor-dependent transcription.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of DMAP1 expression in HeLa cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



产品名称: DMAP1 Rabbit Monoclonal Antibody
产品货号: AMRe84627



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .