

产品名称: NAT10 Rabbit Monoclonal Antibody
产品货号: AMRe84631



产品概述 (Summary)

产品名称 (Production Name)	NAT10 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NAT10
别名 (Alternative Names)	ALP; hALP; N acetyltransferase like protein; NAT10; NET43;;NAT10
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9H0A0.A synthesized peptide derived from human NAT10

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	116 kDa

研究背景 (Background)


产品名称: NAT10 Rabbit Monoclonal Antibody
产品货号: AMRe84631



RNA cytidine acetyltransferase that catalyzes the formation of N4-acetylcytidine (ac4C) modification on mRNAs, 18S rRNA and tRNAs. Catalyzes ac4C modification of a broad range of mRNAs, enhancing mRNA stability and translation. mRNA ac4C modification is frequently present within wobble cytidine sites and promotes translation efficiency.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of NAT10 expression in HeLa cell lysate.

注意事项 (Note)

For research use only .