

产品名称: HTF9C Rabbit Monoclonal Antibody
产品货号: AMRe84638



产品概述 (Summary)

产品名称 (Production Name)	HTF9C Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HTF9C
别名 (Alternative Names)	HTF9C; TRMT2A;;TRMT2A
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q8IZ69.A synthesized peptide derived from human TRMT2A

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 69 kDa ; Observed MW: 75 kDa

研究背景 (Background)

S-adenosyl-L-methionine-dependent methyltransferase that catalyzes the formation of 5-methyl-uridine in tRNAs and


产品名称: HTF9C Rabbit Monoclonal Antibody
产品货号: AMRe84638




some mRNAs. Mainly catalyzes the methylation of uridine at position 54 (m5U54) in cytosolic tRNAs.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of HTF9C expression in A431 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .