

产品名称: QK1 Rabbit Monoclonal Antibody  
产品货号: AMRe84640



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | QK1 Rabbit Monoclonal Antibody         |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB, ICC, IP                            |
| 种属反应性 (Reactivity)     | Human, Mouse                           |

## 产品性能 (Performance)

|                     |                                                                                              |
|---------------------|----------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                 |
| 修饰 (Modification)   | Unmodified                                                                                   |
| 同种型 (Isotype)       | IgG                                                                                          |
| 克隆 (Clonality)      | Monoclonal                                                                                   |
| 形式 (Form)           | Liquid                                                                                       |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.     |
| 储存溶液 (Buffer)       | Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol. |
| 纯化方式 (Purification) | Affinity Purification                                                                        |

## 免疫原信息 (Immunogen)

|                        |                                                     |
|------------------------|-----------------------------------------------------|
| 基因名 (Gene Name)        | QK1                                                 |
| 别名 (Alternative Names) | HKQ; Hqk; HQK1; Hqkl; QK1; QK3; QKI; QKI1;; QKI     |
| 基因 ID (Gene ID)        |                                                     |
| 蛋白 ID (SwissProt ID)   | Q96PU8.A synthesized peptide derived from human QKI |

## 产品应用 (Application)

|                          |                                                |
|--------------------------|------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:1000-1:2000, ICC 1:50-1:200, IP 1:20-1:50 |
| 蛋白分子量 (Molecular Weight) | 38 kDa                                         |

## 研究背景 (Background)

RNA reader protein, which recognizes and binds specific RNAs, thereby regulating RNA metabolic processes, such as pre-


**产品名称: QK1 Rabbit Monoclonal Antibody**  
**产品货号: AMRe84640**




mRNA splicing, circular RNA (circRNA) formation, mRNA export, mRNA stability and/or translation. Involved in various cellular processes, such as mRNA storage into stress granules, apoptosis, lipid deposition, interferon response, glial cell fate and development.

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of QK1 expression in K562 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

## 注意事项 (Note)

For research use only .