

产品名称: Themis Rabbit Monoclonal Antibody
产品货号: AMRe84642



产品概述 (Summary)

产品名称 (Production Name)	Themis Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Themis
别名 (Alternative Names)	Protein THEMIS; SPOT; THEMIS; TSEPA;;THEMIS
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q8N1K5.A synthesized peptide derived from human THEMIS

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	73 kDa

研究背景 (Background)

Plays a central role in late thymocyte development by controlling both positive and negative T-cell selection. Required to


产品名称: Themis Rabbit Monoclonal Antibody
产品货号: AMRe84642



sustain and/or integrate signals required for proper lineage commitment and maturation of T-cells. Regulates T-cell development through T-cell antigen receptor (TCR) signaling and in particular through the regulation of calcium influx and phosphorylation of Erk.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Themis expression in Jurkat cell lysate.

注意事项 (Note)

For research use only .