

产品名称: NCS1 Rabbit Monoclonal Antibody
产品货号: AMRe84658



产品概述 (Summary)

产品名称 (Production Name)	NCS1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NCS1
别名 (Alternative Names)	FLUP; FREQ; Frequentin; Mfreq; ncs1;;NCS 1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P62166.A synthesized peptide derived from human NCS 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	22 kDa

研究背景 (Background)

Neuronal calcium sensor, regulator of G protein-coupled receptor phosphorylation in a calcium dependent manner.


产品名称: NCS1 Rabbit Monoclonal Antibody
产品货号: AMRe84658



Directly regulates GRK1 (RHOK), but not GRK2 to GRK5. Can substitute for calmodulin.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of NCS1 expression in 293T cell lysate.

注意事项 (Note)

For research use only .