

产品名称: Olig3 Rabbit Monoclonal Antibody
产品货号: AMRe84659



产品概述 (Summary)

产品名称 (Production Name)	Olig3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Olig3
别名 (Alternative Names)	bHLHB7; bHLHe20; Olig3; Oligo3;;Oligo3
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q7RTU3.A synthesized peptide derived from human Oligo3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	29 kDa

研究背景 (Background)

May determine the distinct specification program of class A neurons in the dorsal part of the spinal cord and suppress


产品名称: Olig3 Rabbit Monoclonal Antibody
产品货号: AMRe84659



specification of class B neurons.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Olig3 expression in 293T cell lysate.

注意事项 (Note)

For research use only .