

产品名称: USP24 Rabbit Monoclonal Antibody
产品货号: AMRe84661



产品概述 (Summary)

产品名称 (Production Name)	USP24 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	USP24
别名 (Alternative Names)	Usp24;;USP24
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9UPU5.A synthesized peptide derived from human USP24

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	294 kDa

研究背景 (Background)

Ubiquitin-specific protease that regulates cell survival in various contexts through modulating the protein stability of some


产品名称: USP24 Rabbit Monoclonal Antibody
产品货号: AMRe84661



of its substrates including DDB2, MCL1 or TP53. Plays a positive role on ferritinophagy where ferritin is degraded in lysosomes and releases free iron.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of USP24 expression in HeLa cell lysate.

注意事项 (Note)

For research use only .