

产品名称: MED4 Rabbit Monoclonal Antibody
产品货号: AMRe84662



产品概述 (Summary)

产品名称 (Production Name)	MED4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MED4
别名 (Alternative Names)	ARC36; DRIP36; HSPC126; med4; TRAP36; VDRIP;;MED4
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9NPJ6.A synthesized peptide derived from human MED4

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 30 kDa ; Observed MW: 35 kDa

研究背景 (Background)

Component of the Mediator complex, a coactivator involved in the regulated transcription of nearly all RNA polymerase II-


产品名称: MED4 Rabbit Monoclonal Antibody
产品货号: AMRe84662



dependent genes. Mediator functions as a bridge to convey information from gene-specific regulatory proteins to the basal RNA polymerase II transcription machinery. Mediator is recruited to promoters by direct interactions with regulatory proteins and serves as a scaffold for the assembly of a functional preinitiation complex with RNA polymerase II and the general transcription factors.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of MED4 expression in HepG2 cell lysate.

注意事项 (Note)

For research use only .