

产品名称: UPP1 Rabbit Monoclonal Antibody
产品货号: AMRe84664



产品概述 (Summary)

产品名称 (Production Name)	UPP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	UPP1
别名 (Alternative Names)	UDRPASE; UPase 1; UPASE; UPP; Upp1; UrdPase 1;;Uridine phosphorylase 1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q16831.A synthesized peptide derived from human Uridine phosphorylase 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 34 kDa ; Observed MW: 28 kDa

研究背景 (Background)

Catalyzes the reversible phosphorylytic cleavage of uridine to uracil and ribose-1-phosphate which can then be utilized as


产品名称: UPP1 Rabbit Monoclonal Antibody
产品货号: AMRe84664




carbon and energy sources or in the rescue of pyrimidine bases for nucleotide synthesis.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of UPP1 expression in (1) HeLa cell lysate; (2) 293 cell lysate.

产品名称: UPP1 Rabbit Monoclonal Antibody
产品货号: AMRe84664



注意事项 (Note)

For research use only .