

产品名称: UBE4A Rabbit Monoclonal Antibody
产品货号: AMRe84667



产品概述 (Summary)

产品名称 (Production Name)	UBE4A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	UBE4A
别名 (Alternative Names)	UBE4A; UBOX2; UFD2; UFD2b;;UBE4A
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q14139.A synthesized peptide derived from human UBE4A

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100
蛋白分子量 (Molecular Weight)	123 kDa

研究背景 (Background)

Ubiquitin-protein ligase that probably functions as an E3 ligase in conjunction with specific E1 and E2 ligases. May also


产品名称: UBE4A Rabbit Monoclonal Antibody
产品货号: AMRe84667




function as an E4 ligase mediating the assembly of polyubiquitin chains on substrates ubiquitinated by another E3 ubiquitin ligase. Mediates 'Lys-48'-linked polyubiquitination of substrates.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of UBE4A expression in 293 cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: UBE4A Rabbit Monoclonal Antibody
产品货号: AMRe84667



注意事项 (Note)

For research use only .