

产品名称: GLU2B Rabbit Monoclonal Antibody
产品货号: AMRe84677



产品概述 (Summary)

产品名称 (Production Name)	GLU2B Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GLU2B
别名 (Alternative Names)	G19P1; PCLD; PKCSH; PLD1; PRKCSH;;Glucosidase 2 beta
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P14314.A synthesized peptide derived from human Glucosidase 2 beta

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 59 kDa ; Observed MW: 80 kDa

研究背景 (Background)

Regulatory subunit of glucosidase II that cleaves sequentially the 2 innermost alpha-1,3-linked glucose residues from the


产品名称: GLU2B Rabbit Monoclonal Antibody
产品货号: AMRe84677



Glc2Man9GlcNAc2 oligosaccharide precursor of immature glycoproteins. Required for efficient PKD1/Polycystin-1 biogenesis and trafficking to the plasma membrane of the primary cilia.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of GLU2B expression in A431 cell lysate.

注意事项 (Note)

For research use only .