

产品名称: TUG Rabbit Monoclonal Antibody
产品货号: AMRe84690



产品概述 (Summary)

产品名称 (Production Name)	TUG Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TUG
别名 (Alternative Names)	ASPC; ASPCR1; ASPL; ASPS; Aspscr1; RCC17; TUG; UBXD9; UBXN9;;TUG
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9BZE9.A synthesized peptide derived from human TUG

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 60 kDa ; Observed MW: 70-80 kDa

研究背景 (Background)

Tethering protein that sequesters GLUT4-containing vesicles in the cytoplasm in the absence of insulin. Modulates the


产品名称: TUG Rabbit Monoclonal Antibody
产品货号: AMRe84690



amount of GLUT4 that is available at the cell surface (By similarity).Enhances VCP methylation catalyzed by VCPKMT.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of TUG expression in K562 cell lysate.

注意事项 (Note)

For research use only .