

产品名称: CPSF73 Rabbit Monoclonal Antibody
产品货号: AMRe84691



产品概述 (Summary)

产品名称 (Production Name)	CPSF73 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CPSF73
别名 (Alternative Names)	CPSF; cpsf3; CPSF73; YSH1;; CPSF73
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9UKF6.A synthesized peptide derived from human CPSF73

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100
蛋白分子量 (Molecular Weight)	77 kDa

研究背景 (Background)

Component of the cleavage and polyadenylation specificity factor (CPSF) complex that plays a key role in pre-mRNA 3'-end


产品名称: CPSF73 Rabbit Monoclonal Antibody
产品货号: AMRe84691




formation, recognizing the AAUAAA signal sequence and interacting with poly(A) polymerase and other factors to bring about cleavage and poly(A) addition. Has endonuclease activity, and functions as an mRNA 3'-end-processing endonuclease.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of CPSF73 expression in (1) HeLa cell lysate; (2) NIH/3T3 cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: CPSF73 Rabbit Monoclonal Antibody
产品货号: AMRe84691



注意事项 (Note)

For research use only .