

产品名称: SRP19 Rabbit Monoclonal Antibody
产品货号: AMRe84693



产品概述 (Summary)

产品名称 (Production Name)	SRP19 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SRP19
别名 (Alternative Names)	SRP19;;SRP19
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P09132.A synthesized peptide derived from human SRP19

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 16 kDa ; Observed MW: 20 kDa

研究背景 (Background)

Component of the signal recognition particle (SRP) complex, a ribonucleoprotein complex that mediates the cotranslational


产品名称: SRP19 Rabbit Monoclonal Antibody
产品货号: AMRe84693



targeting of secretory and membrane proteins to the endoplasmic reticulum (ER). Binds directly to 7SL RNA. Mediates binding of SRP54 to the SRP complex.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of SRP19 expression in HeLa cell lysate.

注意事项 (Note)

For research use only .