

产品名称: ORP1 Rabbit Monoclonal Antibody
产品货号: AMRe84696



产品概述 (Summary)

产品名称 (Production Name)	ORP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ORP1
别名 (Alternative Names)	ORP1; OSBP8; OSBPL1; OSBPL1A; OSBPL1B;; ORP1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9BXW6.A synthesized peptide derived from human ORP1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 108 kDa ; Observed MW: 50, 108 kDa

研究背景 (Background)

Binds phospholipids; exhibits strong binding to phosphatidic acid and weak binding to phosphatidylinositol 3-


产品名称: ORP1 Rabbit Monoclonal Antibody
产品货号: AMRe84696




phosphate. Stabilizes GTP-bound RAB7A on late endosomes/lysosomes and alters functional properties of late endocytic compartments via its interaction with RAB7A.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of ORP1 expression in 293 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .