

产品名称: CNTN4 Rabbit Monoclonal Antibody
产品货号: AMRe84697



产品概述 (Summary)

产品名称 (Production Name)	CNTN4 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CNTN4
别名 (Alternative Names)	AXCAM; BIG 2; Cntn4; CNTN4A; contactin 4; SCA16;;Contactin 4
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q8I WV2.A synthesized peptide derived from human Contactin 4

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 113 kDa ; Observed MW: 150 kDa

研究背景 (Background)

Contactins mediate cell surface interactions during nervous system development. Has some neurite outgrowth-promoting


产品名称: CNTN4 Rabbit Monoclonal Antibody
产品货号: AMRe84697



activity. May be involved in synaptogenesis.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of CNTN4 expression in Human cerebellum cell lysate.

注意事项 (Note)

For research use only .