

产品名称: SOCS7 Rabbit Monoclonal Antibody
产品货号: AMRe84701



产品概述 (Summary)

产品名称 (Production Name)	SOCS7 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SOCS7 Ash and phospholipase C gamma-binding protein; NAP4; Nck; Nck Ash and phospholipase C binding protein; Nck associated protein 4; NCKAP4; SH2 domain containing SOCS box protein; SOCS4; SOCS6; Socs7; Suppressor of cytokine signaling 7;;SOCS 7
别名 (Alternative Names)	
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O14512.A synthesized peptide derived from human SOCS 7

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	63 kDa

产品名称: SOCS7 Rabbit Monoclonal Antibody
产品货号: AMRe84701




研究背景 (Background)


Regulates signaling cascades probably through protein ubiquitination and/or sequestration. Functions in insulin signaling and glucose homeostasis through IRS1 ubiquitination and subsequent proteasomal degradation.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of SOCS7 expression in HeLa cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .