

产品名称: HOOK2 Rabbit Monoclonal Antibody
产品货号: AMRe84704



产品概述 (Summary)

产品名称 (Production Name)	HOOK2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HOOK2
别名 (Alternative Names)	hHK2; hhook2; HK2; HOOK 2;;HOOK2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q96ED9.A synthesized peptide derived from human HOOK2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 83 kDa ; Observed MW: 85 kDa

研究背景 (Background)

Component of the FTS/Hook/FHIP complex (FHF complex). The FHF complex may function to promote vesicle trafficking


产品名称: HOOK2 Rabbit Monoclonal Antibody
产品货号: AMRe84704



and/or fusion via the homotypic vesicular protein sorting complex (the HOPS complex). Contributes to the establishment and maintenance of centrosome function. May function in the positioning or formation of aggresomes, which are pericentriolar accumulations of misfolded proteins, proteasomes and chaperones. FHF complex promotes the distribution of AP-4 complex to the perinuclear area of the cell.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of HOOK2 expression in HeLa cell lysate.

注意事项 (Note)

For research use only .