

产品名称: NOP10 Rabbit Monoclonal Antibody
产品货号: AMRe84707



产品概述 (Summary)

产品名称 (Production Name)	NOP10 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NOP10
别名 (Alternative Names)	DKCB1; NOLA3; NOP10; NOP10P;;NOP10
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9NPE3.A synthesized peptide derived from human NOP10

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 8 kDa ; Observed MW: 10 kDa

研究背景 (Background)

Required for ribosome biogenesis and telomere maintenance. Part of the H/ACA small nucleolar ribonucleoprotein (H/ACA


产品名称: NOP10 Rabbit Monoclonal Antibody
产品货号: AMRe84707



snoRNP) complex, which catalyzes pseudouridylation of rRNA

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of NOP10 expression in Jurkat cell lysate.

注意事项 (Note)

For research use only .