

产品名称: NPY5R Rabbit Monoclonal Antibody
产品货号: AMRe84708



产品概述 (Summary)

产品名称 (Production Name)	NPY5R Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NPY5R
别名 (Alternative Names)	NPY 5R; NPY R5; NPY Y5; NPY Y5 receptor; NPY5R; NPYR5; NPY5Y5; Y5 receptor;; NPY5R
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q15761.A synthesized peptide derived from human NPY5R

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100
蛋白分子量 (Molecular Weight)	51 kDa

研究背景 (Background)


产品名称: NPY5R Rabbit Monoclonal Antibody
产品货号: AMRe84708



Receptor for neuropeptide Y and peptide YY. The activity of this receptor is mediated by G proteins that inhibit adenylate cyclase activity. Seems to be associated with food intake. Could be involved in feeding disorders.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of NPY5R expression in Ramos cell lysate.

注意事项 (Note)

For research use only .