

产品名称: Gamma glutamyl hydrolase Rabbit Monoclonal Antibody
产品货号: AMRe84712



产品概述 (Summary)

产品名称 (Production Name)	Gamma glutamyl hydrolase Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Gamma glutamyl hydrolase
别名 (Alternative Names)	folylpolyglamaglutamyl hydrolase; gamma glutamyl hydrolase; GGH; ;gamma Glutamyl hydrolase
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q92820.A synthesized peptide derived from human gamma Glutamyl hydrolase

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	36 kDa


产品名称: Gamma glutamyl hydrolase Rabbit Monoclonal Antibody
产品货号: AMRe84712

研究背景 (Background)


Hydrolyzes the polyglutamate sidechains of pteroylpolyglutamates. Progressively removes gamma-glutamyl residues from pteroylpoly-gamma-glutamate to yield pteroyl-alpha-glutamate (folic acid) and free glutamate. May play an important role in the bioavailability of dietary pteroylpolyglutamates and in the metabolism of pteroylpolyglutamates and antifolates.

研究领域 (Research Area)

图片 (Image Data)




All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: Gamma glutamyl hydrolase Rabbit Monoclonal
Antibody
产品货号: AMRe84712



Western blot analysis of Gamma glutamyl hydrolase expression in MCF7 cell lysate.

注意事项 (Note)

For research use only .