

产品名称: TAF15 Rabbit Monoclonal Antibody
产品货号: AMRe84717



产品概述 (Summary)

产品名称 (Production Name)	TAF15 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TAF15
别名 (Alternative Names)	Npl3; RBP56; RNA binding protein 56; TAF; TAF15; TAF2N; TAFII68;;TAF15
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q92804.A synthesized peptide derived from human TAF15

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 62 kDa ; Observed MW: 75 kDa

研究背景 (Background)


产品名称: TAF15 Rabbit Monoclonal Antibody
产品货号: AMRe84717



RNA and ssDNA-binding protein that may play specific roles during transcription initiation at distinct promoters. Can enter the preinitiation complex together with the RNA polymerase II (Pol II).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of TAF15 expression in HeLa cell lysate.

注意事项 (Note)

For research use only .