

产品名称: TOE1 Rabbit Monoclonal Antibody
产品货号: AMRe84731



产品概述 (Summary)

产品名称 (Production Name)	TOE1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TOE1
别名 (Alternative Names)	Toe1;;Target of EGR1 protein 1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q96GM8.A synthesized peptide derived from human Target of EGR1 protein 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 57 kDa ; Observed MW: 60 kDa

研究背景 (Background)

Inhibits cell growth rate and cell cycle. Induces CDKN1A expression as well as TGF-beta expression. Mediates the inhibitory


产品名称: TOE1 Rabbit Monoclonal Antibody
产品货号: AMRe84731



growth effect of EGR1. Involved in the maturation of snRNAs and snRNA 3'-tail processing.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of TOE1 expression in 293T cell lysate.

注意事项 (Note)

For research use only .