

产品名称: p60 katanin Rabbit Monoclonal Antibody
产品货号: AMRe84741



产品概述 (Summary)

产品名称 (Production Name)	p60 katanin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	p60 katanin
别名 (Alternative Names)	Katanin p60 ATPase containing subunit A1; Katanin p60 subunit A1; KATNA1; p60 katanin;; p60 Katanin
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O75449.A synthesized peptide derived from human p60 Katanin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	56 kDa

研究背景 (Background)


产品名称: p60 katanin Rabbit Monoclonal Antibody
产品货号: AMRe84741




Catalytic subunit of a complex which severs microtubules in an ATP-dependent manner. Microtubule severing may promote rapid reorganization of cellular microtubule arrays and the release of microtubules from the centrosome following nucleation.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of p60 katanin expression in HeLa cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .