

产品名称: GRSF1 Rabbit Monoclonal Antibody
产品货号: AMRe84742



产品概述 (Summary)

产品名称 (Production Name)	GRSF1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GRSF1
别名 (Alternative Names)	G rich RNA sequence binding factor 1; G rich sequence binding factor 1; GRSF1;; GRSF 1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q12849.A synthesized peptide derived from human GRSF 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	53 kDa

研究背景 (Background)


产品名称: GRSF1 Rabbit Monoclonal Antibody
产品货号: AMRe84742




Regulator of post-transcriptional mitochondrial gene expression, required for assembly of the mitochondrial ribosome and for recruitment of mRNA and lncRNA.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of GRSF1 expression in 293 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .