

产品名称: ZHX2 Rabbit Monoclonal Antibody
产品货号: AMRe84743



产品概述 (Summary)

产品名称 (Production Name)	ZHX2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ZHX2 AFP regulator 1; AFR1; Alpha fetoprotein regulator 1; RAF; Regulator of AFP;
别名 (Alternative Names)	Transcription factor ZHX2; ZHX2; Zinc finger and homeodomain protein 2; Zinc fingers and homeoboxes 2; Zinc fingers and homeoboxes protein 2;; ZHX2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9Y6X8.A synthesized peptide derived from human ZHX2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 92 kDa ; Observed MW: 110 kDa

产品名称: ZHX2 Rabbit Monoclonal Antibody
产品货号: AMRe84743




研究背景 (Background)


Acts as a transcriptional repressor. Represses the promoter activity of the CDC25C gene stimulated by NFYA. May play a role in retinal development where it regulates the composition of bipolar cell populations, by promoting differentiation of bipolar OFF-type cells.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of ZHX2 expression in HeLa cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .