

产品名称: SLP2 Rabbit Monoclonal Antibody
产品货号: AMRe84746



产品概述 (Summary)

产品名称 (Production Name)	SLP2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ICC,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SLP2 EPB72 like protein 2; HSPC108; Paraprotein target 7; Paratarg 7; SLP2;
别名 (Alternative Names)	Stomatin (EPB72) like 2; Stomatin like 2; Stomatin like protein 2; Stomatin-like protein 2; STOML2;;SLP 2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9UJZ1.A synthesized peptide derived from human SLP 2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	39 kDa

产品名称: SLP2 Rabbit Monoclonal Antibody
产品货号: AMRe84746




研究背景 (Background)


Mitochondrial protein that probably regulates the biogenesis and the activity of mitochondria. Stimulates cardiolipin biosynthesis, binds cardiolipin-enriched membranes where it recruits and stabilizes some proteins including prohibitin and may therefore act in the organization of functional microdomains in mitochondrial membranes.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of SLP2 expression in (1) Jurkat cell lysate; (2) RAW 264.7 cell lysate.



Immunoprecipitate (IP) analysis using the Antibody at 1:50 dilution. (wb at 1:1K dilution)

注意事项 (Note)

For research use only .