

产品名称: Rad51D Rabbit Monoclonal Antibody
产品货号: AMRe84747



产品概述 (Summary)

产品名称 (Production Name)	Rad51D Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Rad51D
别名 (Alternative Names)	BROVCA4; R51H3; Rad51I3; TRAD;;RAD51D
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O75771.A synthesized peptide derived from human RAD51D

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	35 kDa

研究背景 (Background)

Involved in the homologous recombination repair (HRR) pathway of double-stranded DNA breaks arising during DNA


产品名称: Rad51D Rabbit Monoclonal Antibody
产品货号: AMRe84747



replication or induced by DNA-damaging agents. The BCDX2 complex binds single-stranded DNA, single-stranded gaps in duplex DNA and specifically to nicks in duplex DNA.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Rad51D expression in 293 cell lysate.

注意事项 (Note)

For research use only .