

产品名称: ICAM2/CD102 Rabbit Monoclonal Antibody
产品货号: AMRe84749



产品概述 (Summary)

产品名称 (Production Name)	ICAM2/CD102 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ICAM2/CD102
别名 (Alternative Names)	CD102; CD102 antigen; ICAM 2; Intercellular adhesion molecule 2; Ly60;;ICAM 2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P35330(mouse)(human:P13598).A synthesized peptide derived from mouse ICAM 2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 31 kDa ; Observed MW: 60 kDa

产品名称: ICAM2/CD102 Rabbit Monoclonal Antibody
产品货号: AMRe84749




研究背景 (Background)


May play a role in lymphocyte recirculation by blocking LFA1 dependent cell adhesion. Mediates adhesive interactions important for antigen specific immune response, NK cell mediated clearance, lymphocyte recirculation, and other cellular interactions important for immune response and surveillance.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of ICAM2 / CD102 expression in mouse heart cell lysate.



Immunoprecipitate (IP) analysis using the Antibody at 1:50 dilution. (wb at 1:3K dilution)

注意事项 (Note)

For research use only .