

产品名称: Ube2L3/UBCH7 Rabbit Monoclonal Antibody
产品货号: AMRe84755



产品概述 (Summary)

产品名称 (Production Name)	Ube2L3/UBCH7 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Ube2L3/UBCH7
别名 (Alternative Names)	UBCE7; Ubch7; UbcM4; Ube2I3;;UBE2L3
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P68036.A synthesized peptide derived from human UBE2L3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 18 kDa ; Observed MW: 17 kDa

研究背景 (Background)

Ubiquitin-conjugating enzyme E2 that specifically acts with HECT-type and RBR family E3 ubiquitin-protein ligases. Does


产品名称: Ube2L3/UBCH7 Rabbit Monoclonal Antibody
产品货号: AMRe84755



not function with most RING-containing E3 ubiquitin-protein ligases because it lacks intrinsic E3-independent reactivity with lysine: in contrast, it has activity with the RBR family E3 enzymes, such as PRKN, RNF31 and ARIH1, that function like RING-HECT hybrids. Accepts ubiquitin from the E1 complex and catalyzes its covalent attachment to other proteins. In vitro catalyzes 'Lys-11'-linked polyubiquitination.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Ube2L3 / UBCH7 expression in Jurkat cell lysate.

注意事项 (Note)

For research use only .