

产品名称: DENN Rabbit Monoclonal Antibody
产品货号: AMRe84766



产品概述 (Summary)

产品名称 (Production Name)	DENN Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DENN
别名 (Alternative Names)	IG20; MADD; RAB3GEP;;MADD
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q8WXG6.A synthesized peptide derived from human MADD

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 183 kDa ; Observed MW: 230 kDa

研究背景 (Background)

Guanylnucleotide exchange factor that regulates small GTPases of the Rab family. Converts GDP-bound inactive form of


产品名称: DENN Rabbit Monoclonal Antibody
产品货号: AMRe84766




RAB27A and RAB27B to the GTP-bound active forms.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of DENN expression in HeLa cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .