

产品名称: RASSF2 Rabbit Monoclonal Antibody
产品货号: AMRe84771



产品概述 (Summary)

产品名称 (Production Name)	RASSF2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RASSF2
别名 (Alternative Names)	CENP 34; RASFADIN; RASSF2;;RASSF2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P50749.A synthesized peptide derived from human RASSF2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 38 kDa ; Observed MW: 37 kDa

研究背景 (Background)

Potential tumor suppressor. Acts as a KRAS-specific effector protein. May promote apoptosis and cell cycle arrest. Stabilizes


产品名称: RASSF2 Rabbit Monoclonal Antibody
产品货号: AMRe84771



STK3/MST2 by protecting it from proteasomal degradation.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of RASSF2 expression in Jurkat cell lysate.

注意事项 (Note)

For research use only .