

产品名称: N4BP1 Rabbit Monoclonal Antibody
产品货号: AMRe84774



产品概述 (Summary)

产品名称 (Production Name)	N4BP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	N4BP1
别名 (Alternative Names)	N4BP1;;N4BP1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q6A037(mouse)(human:O75113).A synthesized peptide derived from mouse N4BP1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 99 kDa ; Observed MW: 110 kDa

研究背景 (Background)


产品名称: N4BP1 Rabbit Monoclonal Antibody
产品货号: AMRe84774



Potent suppressor of cytokine production that acts as a regulator of innate immune signaling and inflammation. Acts as a key negative regulator of select cytokine and chemokine responses elicited by TRIF-independent Toll-like receptors (TLRs), thereby limiting inflammatory cytokine responses to minor insults.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of N4BP1 expression in PC12 cell lysate.

注意事项 (Note)

For research use only .