

产品名称: RBM3 Rabbit Monoclonal Antibody
产品货号: AMRe84775



产品概述 (Summary)

产品名称 (Production Name)	RBM3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RBM3
别名 (Alternative Names)	Rbm3; RNPL;;RBM3
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P98179.A synthesized peptide derived from human RBM3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	17 kDa

研究背景 (Background)


产品名称: RBM3 Rabbit Monoclonal Antibody
产品货号: AMRe84775



Cold-inducible mRNA binding protein that enhances global protein synthesis at both physiological and mild hypothermic temperatures. Reduces the relative abundance of microRNAs, when overexpressed. Enhances phosphorylation of translation initiation factors and active polysome formation.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of RBM3 expression in Jurkat cell lysate.

注意事项 (Note)

For research use only .