

产品名称: PEF1 Rabbit Monoclonal Antibody
产品货号: AMRe84778



产品概述 (Summary)

产品名称 (Production Name)	PEF1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PEF1
别名 (Alternative Names)	ABP32; PEF; pef1; PEF1A; Peflin;;Peflin
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9UBV8.A synthesized peptide derived from human Peflin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	30 kDa

研究背景 (Background)

Calcium-binding protein that acts as an adapter that bridges unrelated proteins or stabilizes weak protein-protein


产品名称: PEF1 Rabbit Monoclonal Antibody
产品货号: AMRe84778



complexes in response to calcium. Together with PDCD6, acts as a calcium-dependent adapter for the BCR(KLHL12) complex, a complex involved in endoplasmic reticulum (ER)-Golgi transport by regulating the size of COPII coats.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of PEF1 expression in HeLa cell lysate.

注意事项 (Note)

For research use only .