

产品名称: SPINK5/LEKTI Rabbit Monoclonal Antibody
产品货号: AMRe84780



产品概述 (Summary)

产品名称 (Production Name)	SPINK5/LEKTI Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SPINK5/LEKTI
别名 (Alternative Names)	LEKTI; NETS; SPINK5; VAKTI;;SPINK5
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9NQ38.A synthesized peptide derived from human SPINK5

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	121 kDa

研究背景 (Background)

Serine protease inhibitor, probably important for the anti-inflammatory and/or antimicrobial protection of mucous


产品名称: SPINK5/LEKTI Rabbit Monoclonal Antibody
产品货号: AMRe84780




epithelia. Contribute to the integrity and protective barrier function of the skin by regulating the activity of defense-activating and desquamation-involved proteases. Inhibits KLK5, it's major target, in a pH-dependent manner. Inhibits KLK7, KLK14 CASP14, and trypsin.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of SPINK5 / LEKTI expression in HeLa cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: SPINK5/LEKTI Rabbit Monoclonal Antibody
产品货号: AMRe84780



注意事项 (Note)

For research use only .