

产品名称: NUP50 Rabbit Monoclonal Antibody
产品货号: AMRe84785



产品概述 (Summary)

产品名称 (Production Name)	NUP50 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NUP50
别名 (Alternative Names)	NPAP60; NPAP60L; Nup50;;NUP50
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9UKX7.A synthesized peptide derived from human NUP50

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 50 kDa ; Observed MW: 55 kDa

研究背景 (Background)

Component of the nuclear pore complex that has a direct role in nuclear protein import. Actively displaces NLSs from


产品名称: NUP50 Rabbit Monoclonal Antibody
产品货号: AMRe84785




importin-alpha, and facilitates disassembly of the importin-alpha:beta-cargo complex and importin recycling. Interacts with regulatory proteins of cell cycle progression including CDKN1B.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of NUP50 expression in HeLa cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: NUP50 Rabbit Monoclonal Antibody
产品货号: AMRe84785



注意事项 (Note)

For research use only .