

产品名称: DNAJC12 Rabbit Monoclonal Antibody
产品货号: AMRe84789



产品概述 (Summary)

产品名称 (Production Name)	DNAJC12 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DNAJC12
别名 (Alternative Names)	DNAJC12; HPANBH4; JDP1; mJDP1;;DNAJC12
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9UKB3.A synthesized peptide derived from human DNAJC12

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	23 kDa

研究背景 (Background)

This gene encodes a member of a subclass of the HSP40/DnaJ protein family. Members of this family of proteins are


产品名称: DNAJC12 Rabbit Monoclonal Antibody
产品货号: AMRe84789




associated with complex assembly, protein folding, and export. Two transcript variants encoding distinct isoforms have been identified for this gene.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of DNAJC12 expression in A549 cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: DNAJC12 Rabbit Monoclonal Antibody
产品货号: AMRe84789



注意事项 (Note)

For research use only .