

产品名称: CD32A+CD32B+CD32C Rabbit Monoclonal Antibody
产品货号: AMRe84795



产品概述 (Summary)

产品名称 (Production Name)	CD32A+CD32B+CD32C Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD32A+CD32B+CD32C
别名 (Alternative Names)	CD32; CD32A; CD32B; CD32C; CDw32; FCG2; FcGR; FCGR2; FCGR2A; FCGR2A1; FCGR2B; FCGR2C; Fcgr3; Fcgr3a; FcgrII; FcRII a; FcRII b; FcRII c; FcRIIC; IGFR2;;CD32
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P12318/P31994/P31995.A synthesized peptide derived from human CD32

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 34-36 kDa ; Observed MW: 37-43 kDa

产品名称: CD32A+CD32B+CD32C Rabbit Monoclonal Antibody
产品货号: AMRe84795




研究背景 (Background)

Binds to the region of immunoglobulins gamma. Low affinity receptor. By binding to IgG it initiates cellular responses against pathogens and soluble antigens. Promotes phagocytosis of opsonized antigens.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of CD32A + CD32B + CD32C expression in U937 cell lysate.

注意事项 (Note)

For research use only .