

产品名称: Phospho-Nucleolin(T76) Rabbit Monoclonal Antibody
产品货号: AMRe84796



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Nucleolin(T76) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-Nucleolin(T76)
别名 (Alternative Names)	MS1116; NCL; Nucl; Nucleolin; Protein C23;;p-NCL (T76)
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P19338.A synthesized peptide derived from human NCL around the phosphorylation site of T76

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 77 kDa ; Observed MW: 100 kDa

产品名称: Phospho-Nucleolin(T76) Rabbit Monoclonal Antibody
产品货号: AMRe84796




研究背景 (Background)

Nucleolin is the major nucleolar protein of growing eukaryotic cells. It is found associated with intranucleolar chromatin and pre-ribosomal particles. It induces chromatin decondensation by binding to histone H1. It is thought to play a role in pre-rRNA transcription and ribosome assembly. May play a role in the process of transcriptional elongation. Binds RNA oligonucleotides with 5'-UUAGGG-3' repeats more tightly than the telomeric single-stranded DNA 5'-TTAGGG-3' repeats.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Phospho-Nucleolin (T76) expression in (1) 293T cell lysate; (2) 293T cell treated with Lambda Phosphatase lysate.

注意事项 (Note)

For research use only .