

产品名称: GNAI2 Rabbit Monoclonal Antibody
产品货号: AMRe84798



产品概述 (Summary)

产品名称 (Production Name)	GNAI2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GNAI2
别名 (Alternative Names)	GIP; Gnai2; GNAI2B;;GNAI2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P04899.A synthesized peptide derived from human GNAI2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	40 kDa

研究背景 (Background)

Guanine nucleotide-binding proteins (G proteins) are involved as modulators or transducers in various transmembrane


产品名称: GNAI2 Rabbit Monoclonal Antibody
产品货号: AMRe84798



signaling systems. The G(i) proteins are involved in hormonal regulation of adenylate cyclase: they inhibit the cyclase in response to beta-adrenergic stimuli. May play a role in cell division.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of GNAI2 expression in U-87 MG cell lysate.

注意事项 (Note)

For research use only .