

产品名称: AK2 Rabbit Monoclonal Antibody
产品货号: AMRe84801



产品概述 (Summary)

产品名称 (Production Name)	AK2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AK2
别名 (Alternative Names)	Adenylate kinase 2; ADK2; ak2;;Adenylate kinase 2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P54819.A synthesized peptide derived from human Adenylate kinase 2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	26 kDa

研究背景 (Background)


产品名称: AK2 Rabbit Monoclonal Antibody
产品货号: AMRe84801



Catalyzes the reversible transfer of the terminal phosphate group between ATP and AMP. This small ubiquitous enzyme involved in energy metabolism and nucleotide synthesis that is essential for maintenance and cell growth. Plays a key role in hematopoiesis.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of AK2 expression in HepG2 cell lysate.

注意事项 (Note)

For research use only .