

产品名称: C3orf38 Rabbit Monoclonal Antibody
产品货号: AMRe84802



产品概述 (Summary)

产品名称 (Production Name)	C3orf38 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	C3orf38
别名 (Alternative Names)	C3orf38;;C3orf38
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q5JPI3.A synthesized peptide derived from human C3orf38

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 38 kDa ; Observed MW: 37 kDa

研究背景 (Background)


May be involved in apoptosis regulation.

产品名称: C3orf38 Rabbit Monoclonal Antibody
产品货号: AMRe84802




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of C3orf38 expression in U87-MG cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

产品名称: C3orf38 Rabbit Monoclonal Antibody
产品货号: AMRe84802



For research use only .