

产品名称: MYO19 Rabbit Monoclonal Antibody
产品货号: AMRe84804



产品概述 (Summary)

产品名称 (Production Name)	MYO19 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MYO19
别名 (Alternative Names)	myo19; Myohd1; myosin XIX;; Myosin 19
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q96H55.A synthesized peptide derived from human Myosin 19

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 109 kDa ; Observed MW: 100 kDa

研究背景 (Background)

Actin-based motor molecule with ATPase activity that localizes to the mitochondrion outer membrane. Motor protein that


产品名称: MYO19 Rabbit Monoclonal Antibody
产品货号: AMRe84804



moves towards the plus-end of actin filaments. Required for mitochondrial inheritance during mitosis.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of MYO19 expression in K562 cell lysate.

注意事项 (Note)

For research use only .