

产品名称: AdSS2 Rabbit Monoclonal Antibody
产品货号: AMRe84813



产品概述 (Summary)

产品名称 (Production Name)	AdSS2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AdSS2
别名 (Alternative Names)	ADEH; ADSS; AMPase 2;; AdSS 2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P30520.A synthesized peptide derived from human AdSS 2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 50 kDa ; Observed MW: 45 kDa

研究背景 (Background)

Plays an important role in the de novo pathway and in the salvage pathway of purine nucleotide biosynthesis. Catalyzes the


产品名称: AdSS2 Rabbit Monoclonal Antibody
产品货号: AMRe84813




first committed step in the biosynthesis of AMP from IMP.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of AdSS 2 expression in 293 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

产品名称: AdSS2 Rabbit Monoclonal Antibody
产品货号: AMRe84813



注意事项 (Note)

For research use only .